

Systems. A true 12- lead interpretive ECG. Performing fast testing and managing them easily with you at the point-of-care.

- Cardio Universal EMR Interface™
- Cardio Patient Data Management[™] System
- Built-In Networking
- Unlimited Patient Database
- Legendary quality, excellent customer support and training
- · Reliable, user-friendly, compact and durable devices
- Feature rich at the best price
- Made in USA
- FDA Clearance since 1989





Nasiff CardioCard®

PC BASED RESTING BLUETOOTH® ECG SYSTEM

SPECIFICATIONS

CONNECTION WITH PC

POWER SUPPLY	DC 3v from two AA batteries	2110
POWER CONSUMPTION	< 200mA	
ECG MEASUREMENT RESOLUTIONS	2-4 msec, 1.22uV at 10mm/mV gain	
HEART RATE RANGE	30-300bpm	400
LEADS	Standard 12- lead	Continue (philosophilos
LEAD LEAKAGE	< 10uA	mangagaga kaladadada ing mgagagagagagagagagagagagagagagagagagag
ACQUISITION MODE	Simultaneous interpretive 12- lead acquisition	Lineamente Linking
SAMPLING RATE	12- lead sampling rate 250-1000Hz	and the state of t
A/D RATE	1kHz	AMANANANANANANANANANANANANANANANANANANA
TIME BASE	12.5, 25 and 50 mm/sec	
DYNAMIC RANGE	± 5 mV	- Markakahahahahahahahahahaha
FREQUENCY RESPONSE	0.05 to 150 Hz (higher available)	The control of the co
GAINS	2.5, 5, 10, 10/5, 20, 40 mm/mV	traction of the property of the state of the
INPUT IMPEDANCE	>10 ¹² ohms	to an Taranga Taranga
SYSTEM NOISE	<40 uV RTI	atrautes also alle alle a tracticate automobile and

Bluetooth®

FILTERS

ECG INTERPRETATION CLASSES

High performance baseline filter, baseline correction and other filters to minimize motion artifact

Mls, blocks, enlargements, axis and AFib

CMRR >120 dB

NETWORKING Built-in

MEASUREMENTS

HR with global measurements, complete interval values and ST information

DIMENSIONS 2.75" W x 1.25" D x 4.87" H (6.985 cm x 3.175 cm x 12.364 cm)

WEIGHT 0.26 lb (0.12 kg)

PC REQUIREMENTS

Tablet/Laptop/Desktop PC, Intel Compatible, 512MB RAM, 40GB HD, WINDOWS XP through WINDOWS 7, 8, 10 and MAC's with Boot

Camp, Dual Boot,VMware Fusion or Parallels

WARRANTY 2-year limited warranty, support and training

ORDERING INFORMATION CardioResting BT System — CC-ECG BT

FOR MORE INFORMATION CONTACT:

Nasiff Associates, Inc. 841-1 County Route 37, Central Square, NY 13036 phone: 315.676.2346 • fax: 315.676.4711 • toll-free: 866.NASIFFA (866.627.4332) tech support: 315.676.2346 • email: sales@nasiff.com • web: www.nasiff.com

CPT CODE 93000 (additional codes: 93005; 93010)



